



## Trace elements

Na <sub>2</sub> -EDTA	29.00	g
MnSO <sub>4</sub> x n H <sub>2</sub> O	3.26	g
CoCl <sub>2</sub> x 6 H <sub>2</sub> O	1.80	g
ZnSO <sub>4</sub> x 7 H <sub>2</sub> O	1.00	g
NiSO <sub>4</sub> x 6 H <sub>2</sub> O	0.11	g
CuSO <sub>4</sub> x 5 H <sub>2</sub> O	0.10	g
H <sub>3</sub> BO <sub>3</sub>	0.10	g
AlK(SO <sub>4</sub> ) <sub>2</sub> x 12 H <sub>2</sub> O	0.10	g
Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	0.10	g
Na <sub>2</sub> WO <sub>4</sub> x 2 H <sub>2</sub> O	0.10	g
Na <sub>2</sub> SeO <sub>3</sub>	0.05	g
Distilled water	1000.00	ml

First dissolve EDTA x 2Na and adjust pH to 8.0 with KOH. Then add the minerals.