

Main sol. 388

KH ₂ PO ₄	1.50	g
Na ₂ HPO ₄ x 12 H ₂ O	4.20	g
NH ₄ Cl	0.50	g
MgCl ₂ x 6 H ₂ O	0.38	g
CaCl ₂ x 2 H ₂ O	0.06	g
Fe(NH ₄) ₂ (SO ₄) ₂ x 6 H ₂ O	0.04	g
CoCl ₂ x 6 H ₂ O (0.1% w/v)	2.90	ml
Na ₂ MoO ₄ x 2 H ₂ O (0.1% w/v)	2.40	ml
Na ₂ SeO ₃ x 5 H ₂ O (0.01% w/v)	1.70	ml
MnCl ₂ x 4 H ₂ O (0.1% w/v)	2.00	ml
ZnSO ₄ (0.1% w/v)	2.80	ml
Yeast extract	2.00	g
Polypeptone	2.00	g
Starch (soluble)	5.00	g
Sodium resazurin (0.1% w/v)	0.50	ml
Na ₂ CO ₃	1.00	g
L-Cysteine HCl x H ₂ O	1.00	g
Wolin's vitamin solution (10x)	1.00	ml
Distilled water	1000.00	ml

Dissolve ingredients (except carbonate, cysteine and vitamins), bring medium to the boil, then cool to room temperature under 100% N₂ gas atmosphere. Add and dissolve carbonate and cysteine, adjust pH to 7.2, then dispense under 100% N₂ gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave. Add vitamins from an anoxic stock solution prepared under 100% N₂ gas and sterilized by filtration. Adjust pH of complete medium to 7.2, if necessary.