Microorganisms



503b: FWM MEDIUM (PERCHLORATE)

This recipe contains strain-specific modifications for Thauera pivalivorans DSM 14691 *

Final pH: 7.0 - 7.2 Final volume: 1001 ml

Solution A

Solution A	950.00	ml	
Solution B	30.00	ml	
Solution C	10.00	ml	
Solution D	1.00	ml	
Solution E	10.00	ml	

Sparge solution A with 80% N_2 and 20% CO_2 gas mixture for 30 – 45 min to make it anoxic, distribute under same gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave. Solution B is autoclaved separately under 80% N_2 and 20% CO_2 gas atmosphere. Solutions D and E are prepared under 100% N_2 gas and sterilized by filtration. To complete the medium appropriate amounts of solutions B to E are added to the sterile solution A in the sequence as indicated. Adjust pH of complete medium to 7.0 - 7.2, if necessary.

Jointon A		
KH_2PO_4	0.20	g
NH_4CI	0.25	g
NaCl	1.00	g
$MgCl_2 \times 6 H_2O$	0.40	g
KCI	0.50	g
$CaCl_2 \times 2 H_2O$	0.15	g
Na_2SO_4	0.15	g
Modified Wolin's mineral solution	10.00	ml
NaNO ₃	0.85	g
Distilled water	940.00	ml
Solution B Na ₂ CO ₃	1.50	g
		9

Solution C		
Na acotato	1 26	

30.00

ml

Distilled water

^{*} Omit solution E and acetate. Add 0.85 g/l $NaNO_3$ to solution A and use 1.00 g/l pivalic acid as substrate.

Microorganisms

503b: FWM MEDIUM (PERCHLORATE)



Pivalic acid Distilled water	1.00 10.00	g ml
Distilled water	10.00	
Solution D Wolin's vitamin solution (10x)	1.00	ml
Solution E		
Sodium perchlorate	1.22	g
Distilled water	10.00	ml
Modified Wolin's mineral solution (from	medium 141)	
Nitrilotriacetic acid	1.50	g
$MgSO_4 \times 7 H_2O$	3.00	g
$MnSO_4 \times H_2O$	0.50	g
NaCl	1.00	g
FeSO ₄ x 7 H ₂ O	0.10	g
$CoSO_4 \times 7 H_2O$	0.18	g
CaCl ₂ x 2 H ₂ O	0.10	g
$ZnSO_4 \times 7 H_2O$	0.18	g
CuSO ₄ x 5 H ₂ O	0.01	g
$AIK(SO_4)_2 \times 12 H_2O$	0.02	g
H_3BO_3	0.01	g
$Na_2MoO_4 \times 2 H_2O$	0.01	g
$NiCl_2 \times 6 H_2O$	0.03	g
$Na_2SeO_3 \times 5 H_2O$	0.30	mg
$Na_2WO_4 \times 2 H_2O$	0.40	mg
Distilled water	1000.00	ml

First dissolve nitrilotriacetic acid and adjust pH to 6.5 with KOH, then add minerals. Adjust final to pH 7.0 with KOH.

Wolin's vitamin solution (10x) (from medium 120)

Biotin	20.00	mg
Folic acid	20.00	mg
Pyridoxine hydrochloride	100.00	mg
Thiamine HCI	50.00	mg
Riboflavin	50.00	mg
Nicotinic acid	50.00	mg
Calcium D-(+)-pantothenate	50.00	mg
Vitamin B ₁₂	1.00	mg
p-Aminobenzoic acid	50.00	mg
(DL)-alpha-Lipoic acid	50.00	mg

Microorganisms

503b: FWM MEDIUM (PERCHLORATE)



Distilled water 1000.00 ml