Microorganisms



Main sol. 548

| Beef extract | 1.00 | g |
|--|---------|----|
| Glucose | 10.00 | g |
| N-Z amine (enzymatic digest of casein) | 2.00 | g |
| Yeast extract | 1.00 | g |
| Agar | 20.00 | g |
| Distilled water | 1000.00 | ml |

- 1. Adjust pH to 7.3 with NaOH.
- 2. Modified Bennett's agar can also be used. In it, glucose is replaced by glycerol (10.0 g) and N-Z amine by Bacto casein (2.0 g, Difco).