



Main sol. 619

Methanol	10.00	ml
(NH ₄) ₂ HPO ₄	3.00	g
NaCl	1.00	g
MgSO ₄ x 7 H ₂ O	0.20	g
FeSO ₄ x 7 H ₂ O	10.00	mg
MnSO ₄ x 1 H ₂ O	5.00	mg
Vitamin mixture	5.00	ml
Agar	15.00	g
Tap water	1000.00	ml

Sterilize methanol by filtration, and add to other components sterilized by autoclaving.