



Main sol. 652

NaCl	100.00	g
KCl	5.00	g
MgCl ₂ x 6 H ₂ O	5.00	g
NH ₄ Cl	5.00	g
MgSO ₄ x 7 H ₂ O	5.00	g
Yeast extract	5.00	g
Peptone (Difco)	5.00	g
Trace element solution	5.00	ml
Fe(III) citrate (1%)	5.00	ml
Distilled water	990.00	ml

Adjust pH to 7.5-8.0. Boil the medium for few minutes. Fill in tubes with rubber stoppers under nitrogen gas. Sterilize at 121°C for 15 min.