

## Trace element solution SL-11 (from medium 722)

Na <sub>2</sub> -EDTA x 2 H <sub>2</sub> O	5.20	g
FeCl <sub>2</sub> x 4 H <sub>2</sub> O	1.50	g
ZnCl <sub>2</sub>	70.00	mg
MnCl <sub>2</sub> x 4 H <sub>2</sub> O	100.00	mg
H <sub>3</sub> BO <sub>3</sub>	6.00	mg
CoCl <sub>2</sub> x 6 H <sub>2</sub> O	190.00	mg
CuCl <sub>2</sub> x 2 H <sub>2</sub> O	2.00	mg
NiCl <sub>2</sub> x 6 H <sub>2</sub> O	24.00	mg
Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	36.00	mg
Distilled water	1000.00	ml

Dissolve EDTA in 800 ml distilled water, adjust pH to 7 using 2 N NaOH and add ferrous chloride. After ferrous chloride has dissolved add other compounds. Finally adjust pH to 6.0 and bring volume to 1000 ml.