

Main sol. 1078

Bacto brain heart infusion	0.820	g
PPLO broth (GIBCO)	0.870	g
<i>Hank's BSS</i>	50.000	ml
Distilled water	78.000	ml
Porcine serum (heat-inactivated)	30.000	ml
Yeast extract (25% , autoclaved)	3.320	ml
Phenol red solution (1% , filter-sterilized)	0.133	ml
Distilled water (autoclaved)	0.343	ml

Adjust to pH 7.4.

After autoclaving add: #5-8