



Allen's trace element solution (from medium 88)

MnCl ₂ x 4 H ₂ O	180.00	mg
Na ₂ B ₄ O ₇ x 10 H ₂ O	450.00	mg
ZnSO ₄ x 7 H ₂ O	22.00	mg
CuCl ₂ x 2 H ₂ O	5.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	3.00	mg
VO ₂ SO ₄ x 2 H ₂ O	3.00	mg
CoSO ₄ x 7 H ₂ O	1.00	mg
Distilled water	1000.00	ml

Adjust pH of final solution to 2 with 1 N HCl.