

Main sol. 798

Na ₂ CO ₃	8.00	g
KCl	0.20	g
NaCl	10.00	g
NH ₄ Cl	0.50	g
K ₂ HPO ₄	0.20	g
Yeast extract	4.00	g

Trace element solution SL-10	1.00	ml
-------------------------------------	------	----

Wolin's vitamin solution	1.00	ml
---------------------------------	------	----

Na ₂ S x 9 H ₂ O	0.50	g
--	------	---

Distilled water	1000.00	ml
-----------------	---------	----

1. pH 9.0
2. Prepare the medium anaerobically without the vitamins. Boil the medium, cool under N₂ and dispense under N₂ and autoclave the medium. Add the vitamin solution to the cooled, autoclaved medium.