

Main sol. 828

Potassium phosphate buffer, 1M, pH 7	27.50	ml
MgSO ₄ x 7 H ₂ O	0.05	g
CaCl ₂ x 2 H ₂ O	14.00	mg
Na-lactate	0.56	g
SL-9	0.50	ml
Seven vitamins solution	0.50	ml
Agar, if necessary	15.00	g
Distilled water	972.50	ml

1. Adjust the pH to 7.0 - 7.2.
2. Filter sterilize the vitamin solution and add to the autoclaved medium.