



Main sol. 888

(NH ₄) ₂ SO ₄	1.00	g
MgSO ₄ x 7 H ₂ O	1.00	g
CaCl ₂ x 6 H ₂ O	30.00	mg
Na ₂ HPO ₄	10.00	mg
Sodium succinate	1.00	g
Casamino acids (DIFCO)	1.50	g
Trace elements solution (Pfennig & Lippert,1966)	1.00	ml
Agar (DIFCO)	0.50	g
Distilled water	1000.00	ml

1. pH 7.5

2. Add after autoclaving sterile solutions of:

Na ₂ S ₂ O ₃ x 5 H ₂ O (10% solution)	1.00	ml
Vitamine mixture	1.00	ml