



## Main sol. 911

<b>Bacto Middelbrook 7H11 agar (Difco 283810)</b>	21.00	g
Glycerol	5.00	ml
Distilled water	900.00	ml
<b>Bacto Middelbrook OADC enrichment (Difco 212240)</b>	100.00	ml

1. Bacto Middelbrook OADC enrichment (Difco 212240):
2. Sterilize for 10 min by 121°C cool to 50 - 55°C. Add 100 ml Bacto Middelbrook OADC enrichment Difco (see below) under aseptic conditions and mix thoroughly.
3. Final pH at 25°C 6.6 ± 0.2.