



Trace elements

FeSO ₄ x 7 H ₂ O	2.10	g
Na ₂ -EDTA	5.20	g
H ₃ BO ₃	30.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
ZnSO ₄ x 7 H ₂ O	144.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
Distilled water	1000.00	ml

Adjust the pH to 6.0