



Main sol. 1119

Sodium gluconate	3.00	g
Yeast extract (Difco)	1.00	g
(NH ₄) ₂ SO ₄	1.00	g
KH ₂ PO ₄	1.00	g
NaCl	0.20	g
MgCl ₂ x 6 H ₂ O	0.20	g
CaCl ₂ x 2 H ₂ O	0.05	g
Trace element solution SL8	1.00	ml
Distilled water	1000.00	ml
Agar, for solid medium	1.80	%

1. Adjust pH to 4.0 - 4.5.
2. For solid medium, add 1.8% agar (final) and adjust pH to 5.0 to avoid hydrolysis of agar.