

## Main sol. 1130

NaCl	175.00	g
MgCl <sub>2</sub> x 6 H <sub>2</sub> O	20.00	g
K <sub>2</sub> SO <sub>4</sub>	5.00	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.10	g
Yeast extract (ThermoFischer)	5.00	g
Distilled water	1000.00	ml
Agar (optional)	15.00	g

1. Adjust to pH 6.8-7.2
2. Note: for Agar plates, autoclave NaCl separately. Otherwise the salts will precipitate and the plates not solidify properly.
3. Note: especially for strains DSM 17983 and 17984, yeast extract from ThermoFisher has to be used. Yeast extract from other vendors might not support growth.