

## Main sol. 1133

|                        |         |    |
|------------------------|---------|----|
| Tryptone               | 2.00    | g  |
| Beef extract           | 0.50    | g  |
| Yeast extract          | 0.50    | g  |
| Sodium acetate         | 0.20    | g  |
| Agar                   | 15.00   | g  |
| Distilled water        | 1000.00 | ml |
| Agar, for solid medium | 4.00    | g  |

1. Adjust to pH 7.2-7.4 and autoclave at 121°C/15 min.
2. For soft agar for the maintenance of active cultures reduce agar content to 4.0 g/l.  
Dispense in 3 to 4 ml amounts in 7 ml screw-capped bottles (bijou bottles).