

Wolfe's Mineral Solution

Nitrilotriacetic acid	1.50	g
MgSO ₄	3.00	g
MnSO ₄	0.50	g
NaCl	1.00	g
FeSO ₄	0.10	g
CaCl ₂	0.10	g
CoCl ₂	0.10	g
ZnSO ₄	0.10	g
CuSO ₄	0.01	g
AlK(SO ₄) ₂	0.01	g
H ₃ BO ₃	0.01	g
Na ₂ MoO ₄	0.01	g
Distilled water	1000.00	ml

1. First dissolve nitrilotriacetic acid and adjust pH to 6.5 with KOH, then add minerals.
2. Adjust to pH 8.0 with KOH and sterilize by autoclaving.