



Main sol. 1396

L-Glutamic acid	2.00	g
Trisodium citrate x 2 H ₂ O	2.00	g
KH ₂ PO ₄	1.00	g
NH ₄ Cl	1.00	g
K ₂ SO ₄	5.00	g
NaCl	200.00	g
Casamino acids	4.00	g
Yeast extract	2.00	g
FeSO ₄ x 6 H ₂ O	4.00	mg
Trace Element Solution	2.00	ml
Agar	20.00	g
Distilled water	800.00	ml

Adjust to pH 4.5 with 40% KOH and sterilize by autoclaving.