



Main sol. 1431

Yeast extract	10.00	g
Casamino acids	7.50	g
Trisodium citrate	3.00	g
KCl	2.00	g
MgSO ₄ x 7 H ₂ O	1.00	g
MnCl ₂ x 4 H ₂ O	0.36	mg
FeSO ₄ x 7 H ₂ O	50.00	mg
NaCl	200.00	g
Na ₂ CO ₃	50.00	g
Agar	2.00	%

Liquid medium may be prepared by adding all components together. Agar (2%) medium should be prepared without Na₂CO₃ and sterilised by autoclaving. Add 20-30 ml of a sterile 10% Na₂CO₃ solution to the autoclaved medium containing agar so that the final pH is 9.0.