



Main sol. 1492

| | | |
|--|----------|----|
| Yeast extract | 0.250 | g |
| Proteose peptone (Difco no. 3) | 0.250 | g |
| Casamino acids | 0.250 | g |
| Glucose | 0.250 | g |
| Starch (soluble) | 0.250 | g |
| Na-pyruvate | 0.250 | g |
| K ₂ HPO ₄ | 0.150 | g |
| MgSO ₄ x 7 H ₂ O | 0.025 | g |
| Distilled water | 1000.000 | ml |

1. Adjust pH to 7.2 with K₂HPO₄ or KH₂PO₄.
2. Add agar after adjusting the pH, then boil to dissolve the agar

| | | |
|------|--------|---|
| Agar | 15.000 | g |
|------|--------|---|

3. Autoclave at 121 °C for 15 min.