



## Trace elements III

$\text{NiSO}_4 \times 6 \text{H}_2\text{O}$	13.10	mg
$\text{SnCl}_2 \times 2 \text{H}_2\text{O}$	11.80	mg
$\text{NaVO}_3$	61.00	mg
$(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \times 4 \text{H}_2\text{O}$	1240.00	mg
HCl (10 mM)	100.00	ml

Prepare Trace elements III, flush with nitrogen to make them anoxic, filter-sterilize and store at 4°C