



Mineral Salts Solution for liquid Medium

L-Glutamic acid	1.00	g
Trisodium citrate x 2 H ₂ O	1.00	g
KH ₂ PO ₄	0.50	g
NH ₄ Cl	0.50	g
K ₂ SO ₄	2.50	g
MgCl ₂ x 7 H ₂ O	20.00	g
NaCl	100.00	g
Casamino acids	2.00	g
Yeast extract	1.00	g
FeSO ₄	2.00	mg
Trace Element Solution	1.00	ml
Distilled water	500.00	ml

Adjust to pH 5.5 with 10% KOH and sterilize by autoclaving.