

Main sol. 1536

KH ₂ PO ₄	0.950	g
K ₂ HPO ₄	0.600	g
NH ₄ Cl	0.270	g
Na-propionate	0.480	g
MgSO ₄ x 7 H ₂ O	0.025	g
CaCl ₂ x 2 H ₂ O	0.010	g
Vitamine solution	1.000	ml
Chelated iron solution	1.000	ml
Trace element solution	0.100	ml
Distilled water	1000.000	ml

1. The pH of the final solution is adjusted to 6.3.
2. Note: Please autoclave the chelated iron solution before you add to the medium, so that solves everything!!!!