



## Main sol. 1578

Yeast extract	0.05	g
Proteose peptone (Difco no. 3)	0.05	g
Casamino acids	0.05	g
Glucose	0.05	g
Starch (soluble)	0.05	g
Na-pyruvate	0.03	g
K <sub>2</sub> HPO <sub>4</sub>	0.03	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	5.00	mg
Agar	15.00	g
Distilled water	1000.00	ml

Final pH 7.2; adjust with crystalline K<sub>2</sub>HPO<sub>4</sub> or KH<sub>2</sub>PO<sub>4</sub> before adding agar. Add agar, heat medium to boiling to dissolve agar, and autoclave for 15 min at 121°C.