



Main sol. J358

Trace element solution SLA	1.00	ml
Vitamin solution VA	1.00	ml
KH ₂ PO ₄	1.00	g
MgCl ₂ x 6 H ₂ O	0.50	g
CaCl ₂ x 2 H ₂ O	0.10	g
NH ₄ Cl	1.00	g
NaCl	1.00	g
Sodium acetate	1.00	g
Sodium succinate	1.00	g
Yeast extract (BD-Difco)	0.50	g
Bacto peptone (BD-Difco)	0.30	g
Distilled water	1000.00	ml

1. Adjust pH to 6.9.
2. Separately sterilize vitamin solution VA by filtration, and aseptically add to the medium.