



Main sol. J483

Glucose	10.00	g
Tryptone (BD-Difco)	4.00	g
Beef extract (BD-Difco)	2.00	g
Yeast extract (BD-Difco)	1.00	g
NaCl	4.00	g
K ₂ HPO ₄	1.50	g
FeSO ₄ x 7 H ₂ O	0.10	g
MgCl ₂ x 6 H ₂ O	0.20	g
L-Cysteine HCl x H ₂ O	0.50	g
Microelement solution	10.00	ml
Vitamin solution	10.00	ml
Distilled water	1000.00	ml

Mix components except vitamin solution, and adjust pH to 4.5 - 6.5. After autoclaving, aseptically add filter-sterilized vitamin solution.