



SL-7 trace element solution

HCl (7.7 M, 25%)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
CoCl ₂ x 6 H ₂ O	190.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
ZnCl ₂	70.00	mg
H ₃ BO ₃	62.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg
CuCl ₂ x 2 H ₂ O	17.00	mg
Distilled water	1000.00	ml

Initially, dissolve FeCl₂ x 4H₂O in 25% HCl, then add distilled water followed by the other components.