

Main sol. 205

Casamino acids	15.00	g
Trisodium citrate x 2 H ₂ O	3.00	g
Glutamic acid	2.50	g
MgSO ₄ x 7 H ₂ O	2.50	g
KCl	2.00	g
NaCl	250.00	g
Agar	20.00	g
Distilled water	1000.00	ml

Dissolve agar by heating before adding sodium chloride. Adjust pH to 7.0 before autoclaving. Adjust pH to 8.5 with sterile 5% Na₂CO₃ after heat sterilization of the medium.