Microorganisms



Main sol. 246 (from medium 246)

Beef extract	10.00	g
Peptone	10.00	g
Agar	20.00	g
Tap water	250.00	ml
Sea water	750.00	ml

- 1. Dissolve beef extract and peptone by heating in tap water, adjust pH to 7.8 and boil for 10 min. Readjust pH to 7.3. Add agar and autoclave at 121°C for 20 min. Directly after autoclaving, add warm (55°C) sterile sea water. Liquid medium without agar should be combined when cooled to room temperature.
- 2. Natural sea water is stored in the dark for at least three weeks to "age". If natural sea water is not available use artificial sea water.