



Trace element solution (A) (from medium 253)

FeCl ₂ x 4 H ₂ O	1.80	g
CoCl ₂ x 6 H ₂ O	250.00	mg
NiCl ₂ x 6 H ₂ O	10.00	mg
CuCl ₂ x 2 H ₂ O	10.00	mg
MnCl ₂ x 4 H ₂ O	70.00	mg
ZnCl ₂	100.00	mg
H ₃ BO ₃	500.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	30.00	mg
Na ₂ SeO ₃ x 5 H ₂ O	10.00	mg
Distilled water	1000.00	ml

For dissolving adjust pH to about 3 with 1 N HCl.