



## M 1757

<b>Salt mix 10x</b>	100.00	ml
<b>CasAA mix 5x</b>	200.00	ml
<b>Trace salt mix 100x</b>	10.00	ml
<b>Iron stock solution 2.0 g/L</b>	3.00	ml
<b>Vitamine mixture</b>	5.00	ml
<b>Glucose stock 20%</b>	10.00	ml

1. fill to 1 Liter with H<sub>2</sub>O (autoclaved)
2. anaerobize with nitrogen, can be stored anaerobically