



Mail sol. C₉₂

K ₂ HPO ₄ (1 g/l stock solution)	20.00	ml
MgSO ₄ x 7 H ₂ O (1 g/l stock solution)	20.00	ml
KNO ₃ (10 g/l stock solution)	20.00	ml
SE2 (Soil Extract 2)	100.00	ml
Deionized water	1000.00	ml

1. Make up to 1 litre with deionized water and autoclave at 15 psi for 15 minutes.
2. At the CCAP, SE1 is used for marine algae, SE2 for freshwater and terrestrial protozoa.