

Trace element solution (SLA) (from medium 351)

FeCl ₂ x 4 H ₂ O	1.80	g
CoCl ₂ x 6 H ₂ O	250.00	mg
NiCl ₂ x 6 H ₂ O	10.00	mg
CuCl ₂ x 2 H ₂ O	10.00	mg
MnCl ₂ x 4 H ₂ O	70.00	mg
ZnCl ₂	100.00	mg
H ₃ BO ₃	500.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	30.00	mg
Na ₂ SeO ₃ x 5 H ₂ O	10.00	mg
Distilled water	1000.00	ml

Adjust pH to 2.0 - 3.0.