Microorganisms



104: PYG MEDIUM (modified)

This recipe contains strain-specific modifications for Leptotrichia shahii DSM 19757 *

Final pH: 7.2 Final volume: 1000 ml

Trypticase peptone	5.00	g
Peptone	5.00	g
Yeast extract	10.00	g
Beef extract	5.00	g
Glucose	5.00	g
K ₂ HPO ₄	2.00	g
Tween 80	1.00	ml
Salt solution	40.00	ml
Resazurin	1.00	mg
Vitamin K ₁ solution	0.20	ml
Haemin solution	10.00	ml
L-Cysteine HCl x H_2O	0.50	g
Distilled water	950.00	ml

The vitamin K_1 , haemin solution and the cysteine are added after the medium has been boiled and cooled under CO₂. Adjust pH to 7.2 using 8 N NaOH. Distribute under N₂ and autoclave.

* Anaerobic (90% N₂, 5% H₂, 5% CO₂)

0.25	g
0.50	g
1.00	g
1.00	g
10.00	g
2.00	g
1000.00	ml
	0.50 1.00 1.00 10.00 2.00

Haemin solution (from medium 78)		
Haemin	50.00	mg
NaOH (1 N)	1.00	ml
Distilled water	100.00	ml

Dissolve 50 mg haemin in 1 ml 1 N NaOH; make up to 100 ml with distilled water and filter sterilize. Store refrigerated.

Microorganisms

104: PYG MEDIUM (modified)



Vitamin K₁ solution (from medium 78)

Vitamin K ₁	0.10	ml
Ethanol (95 %)	20.00	ml

Dissolve 0.1 ml of vitamin K_1 in 20 ml 95% ethanol and filter sterilize. Store refrigerated in a brown bottle.