Microorganisms



104b: PY + X MEDIUM

This recipe contains strain-specific modifications for Defluviitalea raffinosedens DSM 28090 *

Final pH: * 8.5

Final volume: 1000 ml

Trypticase peptone (BD BBL)	5.00	g
Meat peptone (pepsin-digested)	5.00	g
Yeast extract	10.00	g
Salt solution	40.00	ml
Sodium resazurin (0.1% w/v)	0.50	ml
L-Cysteine HCl x H ₂ O	0.50	g
D-Glucose	5.00	g
Distilled water	960.00	ml

Dissolve ingredients (except cysteine and glucose) and sparge medium with 100% N_2 gas for 30 - 45 min to make it anoxic. Add cysteine and adjust pH to 7.0, then dispense medium under 100% N_2 gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave. After autoclaving add glucose or any other substrate from sterile anoxic stock solutions prepared under 100% N_2 gas atmosphere. Adjust pH of complete medium to 6.8 - 7.0, if necessary.

Salt solution (from medium 104)

0.25	g
0.50	g
1.00	g
1.00	g
10.00	g
2.00	g
1000.00	ml
	0.50 1.00 1.00 10.00 2.00

^{*} Adjust pH to 8.5 with a sterile anoxic stock solution of Na_2CO_3 (5% w/v) prepared under 80% N_2 and 20% CO_2 gas atmosphere.