

1429: SABOURAUD GLUCOSE MEDIUM

This recipe contains strain-specific modifications for *Acidomonas methanolica* DSM 5432 *

Final pH: * 4.0 - 5.0

Final volume: 1000 ml

SABOURAUD- Glucose-Bouillon (Merck 108339, 2 %)	30.00	g
Agar	15.00	g
Acetic acid	1.00	ml/l
Distilled water	1000.00	ml

* After autoclaving add 1.0 ml/L of filter sterilized acetic acid (99% w/v). The pH should be 4-5.