

1429: SABOURAUD GLUCOSE MEDIUM

This recipe contains strain-specific modifications for *Komagataeibacter europaeus* DSM 6160 *

Final pH: * 3.0 - 4.0

Final volume: 1000 ml

SABOURAUD- Glucose-Bouillon (Merck 108339, 2 %)	30.00	g
Agar	15.00	g
Acetic acid	3.20	ml/l
Distilled water	1000.00	ml

* After autoclaving add 3.2 ml/L of filter sterilized acetic acid (99% w/v). The pH should be 3-4