

1713: GAM Broth (Gifu Anaerobic Medium Broth)

This recipe contains strain-specific modifications for *Mogibacterium kristiansenii* DSM 106282 *

Final volume: 1000 ml

GAM Broth powder (HyServe)	59.00	g
Sodium pyruvate	100.00	mM
Sodium propionate	2.00	M
Sodium acetate	4.50	M
Deionized water	1000.00	ml

Suspend the dehydrated medium in deionized water, and heat to dissolve the medium, sterilize by autoclaving at 115 °C for 15 minutes.

* Add 100 mM sodium pyruvate, 4,5 M sodium acetate and 2 M sodium propionate; anaerobic