

190: YpSs MEDIUM

This recipe contains strain-specific modifications for *Dispira cornuta* DSM 3305 *

Final pH: 6.5

Final volume: 1000 ml

Yeast extract	4.00	g
Starch (soluble)	15.00	g
K ₂ HPO ₄	1.00	g
MgSO ₄ x 7 H ₂ O	0.50	g
Agar	15.00	g
Distilled water	1000.00	ml

pH to 6.5. Sterilize at 121°C for 15 min.

* Coculture with *Cokeromyces recurvatus*