

1: NUTRIENT AGAR

This recipe contains strain-specific modifications for *Paenibacillus agarexedens* DSM 1475 *

Final pH: 7.0

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, for solid medium	15.00	g
Urea	1.00	g/l
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. For Bacillus strains the addition of 10.0 mg $\text{MnSO}_4 \times \text{H}_2\text{O}$ is recommended for sporulation.

* Plus 0.1% urea