

1: NUTRIENT AGAR

This recipe contains strain-specific modifications for *Xanthomonas vesicatoria* DSM 50861 *

Final pH: 7.0

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, for solid medium	15.00	g
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. For *Bacillus* strains the addition of 10.0 mg $\text{MnSO}_4 \times \text{H}_2\text{O}$ is recommended for sporulation.

* Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth medium 1. Subsequent subculturing may be carried out in liquid medium or on agar plates