

1: NUTRIENT AGAR

This recipe contains strain-specific modifications for *Bacillus* sp. DSM 5979 *

Final pH: 7.0

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, for solid medium	15.00	g
Yeast extract	4.00	g/l
Na-benzoate	2.00	g/l
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. For *Bacillus* strains the addition of 10.0 mg $\text{MnSO}_4 \times \text{H}_2\text{O}$ is recommended for sporulation.

* With 0.4% yeast extract and 0.2% sodiumbenzoate