

1: NUTRIENT AGAR

This recipe contains strain-specific modifications for *Bacillus licheniformis* DSM 7459 *

Final pH: 7.0

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, for solid medium	15.00	g
NiCl ₂	200.00	mg/l
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. For *Bacillus* strains the addition of 10.0 mg MnSO₄ × H₂O is recommended for sporulation.

* Plus 200 mg/l NiCl₂