

1: NUTRIENT AGAR

This recipe contains strain-specific modifications for *Methylobacterium mesophilicum* DSM 1708 *

Final pH: 7.0

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, for solid medium	15.00	g
Methanol	1.00	%(v/v)
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. For Bacillus strains the addition of 10.0 mg $\text{MnSO}_4 \times \text{H}_2\text{O}$ is recommended for sporulation.

* With 1% v/v methanol