

## 220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Alteribacter populi* DSM 104632 \*

Final pH: \* 9.0

Final volume: 1000 ml

|                 |         |    |
|-----------------|---------|----|
| Casein peptone  | 15.00   | g  |
| Soy peptone     | 5.00    | g  |
| NaCl            | 5.00    | g  |
| Agar            | 15.00   | g  |
| Distilled water | 1000.00 | ml |

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

\* Adjust pH of final medium to 9.0.