

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Ureibacillus thermosphaericus* DSM 12655 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Soil extract	50.00	%
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus 50% soil extract