

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Alkalicoccus saliphilus* DSM 15402 *

Final pH: * 9.0

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
NaCl	160.00	g/l
MnSO ₄	10.00	mg/l
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus 16% NaCl and 10 mg/l MnSO₄; pH 9