

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Halalkalibacter hemicellulosilyticus* DSM 16731 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus 1/10 Na-sesquicarbonate solution; (see medium 31)