

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Gracilibacillus boracitolerans* DSM 17256 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
MnSO ₄	10.00	mg/l
H ₃ BO ₃	0.05	M
NaCl	2.50	%
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus MnSO₄ (10µg/ml), 0.05M H₃BO₃, 2.5 % NaCl